

Technical Specifications for the Sulphide Ion-Selective Electrode ELIT 8225

Introduction

The Sulphide Ion-Selective Electrode has a solid-state crystal membrane. The electrode is designed for the detection of sulphide ions (S^{2-}) in aqueous solutions and is suitable for use in both field and laboratory applications. The Sulphide Ion is a bivalent anion.

One mole of (S^{2-}) has a mass of 32.064 grams; 1000 ppm is 0.0312M

Dissolve 7.491g sodium sulphide nonahydrate ($Na_2S \cdot 9H_2O$) in 1 litre water.

Physical Specifications

| | |
|----------------------------------|------------|
| Length of body excl gold contact | 130 mm |
| Length of body incl gold contact | 140 mm |
| Diameter of body | 8 mm |
| DC resistance at 25°C | < 2.5 MOhm |
| Minimum feasible sample volume | 5 ml |

Chemical / Operational Specifications

| | |
|--|---|
| Preconditioning / Standard solution | Aprox. 1300 ppm S^{2-} as Na_2S (Best diluted from saturated solution and titrated with Lead Perchlorate standard) |
| Preconditioning time | 5 minutes |
| Optimal pH range | pH 13 to pH 14 |
| Temperature range | 0 to 80 °C |
| Recommended ISAB | 10M NaOH (add 2% v/v) |
| Recommended reference electrode | Double junction (ELIT 003) |
| Reference electrode outer filling solution | 0.1M CH_3COOLi |
| Electrode slope at 25°C | 26 ± 3 mV/decade |
| Concentration range | 0.003 to 3,200 ppm (9×10^{-8} to 0.1 Molar) |
| Response time | < 10 seconds (Defined as time to complete 90% of the change in potential after immersion in the new solution.) |
| Potential drift (in 1000 ppm) | < 3 mV/day (8 hours) (Measured at constant temperature and with ISE and Reference Electrode continually immersed) |

Interference:

Mercury or Silver ions have very high interference and can only be tolerated in very low concentrations relative to the Sulphide – ideally they should be absent.

Note high and narrow pH range (13 to 14).

WARNING: 10M NaOH is a very caustic solution and should be handled with care.

For more information, see: www.nico2000.net.